

EC-12R RESTING ECG SYSTEM (USB/BLUETOOTH)

- SIMPLE, USER FRIENDLY SOFTWARE WITH MULTIPLE FUNCTIONS
- PRECISE QRS CLASSIFICATION
- BUILT-IN EXPERT SYSTEM,
- SHORT AND LONG RESTING ECG RECORD
- VARIOUS PRINTABLE REPORTS,
- CUSTOMIZABLE ALARM CRITERIA
- COMPARING MEASUREMENTS
- DICOM, GDT, MFER, HL7
- LOCAL AND NETWORK OPERATION
- EXPORT AND IMPORT FUNCTIONS
- ARRHYTHMIA ANALYSIS, ARRHYTHMIA OVERVIEW

Monitoring patient with continuous recording


Besides conventional or short-term resting ECG records, you can also make records lasting several hours.

ECG Filtering to improve signal quality


Using an effective digital filtering technology, there is practically no baseline movement. The smoothing filter removes muscle and network noise components from the signals without distorting QRS complexes.

Rhythm view

All resting/continuous/stress test ECG examination will have a window below the normal 1-minute-long ECG display. This full view display shows a 5-minute-long ECG section and we can keep track of all arrhythmias in this section.



MAIN WINDOW



CARDIOSPY LITE